

**EPIPEN® EPIPEN Jr®**

(Adrenalin) Auto-Injektoren 0,3/0,15 mg

## Certificate for Airline Passenger (for attention of airport security and control)

### Statement for rescue medication

**EpiPen®/EpiPen® Junior** is used for the emergency treatment of sudden and life-threatening allergic reactions (anaphylactic shock) to insect stings and bites, foods, drugs, latex and exercise. EpiPen® contains a sterile solution of adrenaline (epinephrine) for emergency intra muscular self-injection.

The patient \_\_\_\_\_ (name)

is suffering from a life threatening \_\_\_\_\_ allergy.

The EpiPen®/EpiPen® Junior (adrenaline) auto-injector is the prescribed rescue medication that the patient has to **carry at all-times**. In case of an emergency (a life threatening anaphylactic shock or a serious allergic reaction) the adrenaline can be administered and further damage caused to the patient's health can be avoided.

Physician  
stamp – signature – date

### Symptome die eine Anaphylaxie ankündigen:

- Kribbeln und Schwellungen an Lippen, Augen, im Rachenraum und im Gesicht
- Hautausschlag / Juckreiz
- Atemnot
- Blutdruckabfall

# EPIPEN® EPIPEN Jr®

(Adrenalin) Auto-Injektoren 0,3/0,15 mg

## Anwendung von EpiPen®/EpiPen® Junior:

1. Umfassen Sie den EpiPen mit der Faust. Der Daumen muss neben dem grauen Sicherheitsdeckel zu liegen kommen.
2. Mit der anderen Hand den grauen Sicherheitsdeckel entfernen.
3. Mit dem schwarzen Ende voraus, den EpiPen in Richtung des äusseren Oberschenkels halten.
4. Den Auto-Injektor im 90° Winkel so kräftig gegen den äusseren Oberschenkel stossen, dass der Injektionsmechanismus ausgelöst wird.
5. Nach Eindringen der Nadel den Auto-Injektor 10 Sekunden in dieser Stellung halten.
6. Auto-Injektor wegnehmen und Einstichstelle 10 Sekunden massieren.
7. Sofort einen Arzt oder eine Ärztin aufsuchen.



## Anwendung eines zweiten EpiPen®/EpiPen® Junior:

Eine zweite Dosis kann, wenn nötig, nach 5–15 Minuten appliziert werden.

**Verfalldatum beachten:** Den EpiPen®/EpiPen® Junior unbedingt vor Erreichen des Verfalldatums (siehe «EXP» auf Packung oder Pen) erneuern. Registrieren Sie sich beim kostenlosen Erinnerungs-Service unter: [www.epipen.ch](http://www.epipen.ch).

**MEDA**

MEDA Pharma GmbH  
Hegnaustrasse 60  
8602 Wangen-Brüttsellen  
Tel. +41 44 835 26 26  
Fax +41 44 835 26 27  
[www.medapharma.ch](http://www.medapharma.ch)

